

BULLET ICE MAKER LTICE8-1



BULLET ICE MAKER LTICE8-1

Used in Beverage industries, Biochemical and biological studies, Refrigerated perishables, Food industries, Medical laboratories, Healthcare, Pharmaceutical.

Also known as Ice Maker, Laboratory Bullet Ice Maker.

LTICE8-1 BULLET ICE MAKER

Vertical stirring engine model

Refrigeration system produces uniform shape ice bullets

Bell-shaped ice with a hexagon hole

Highly durable stainless steel 304 2B built-up

Automatic controller increases working reliability

Large ice storage bin

Eco-friendly, energy saving non-Freon refrigerant



SPECIFICATIONS

Model	LTICE8-1
Capacity	20 kg/24 h
Storage Bin Capacity	6 kg
Refrigerant	R134a
Condenser Type	Air Cooled
Overall Dimension	330x470x605 mm
Packing Dimension	405x555x665 mm
Weight	28/32 kg
Power	220 W
Power Supply	220V /110 V, 50 / 60 Hz

OPTIONAL ACCESSORIES

Accessory Code	Name
2900808008	Water filter



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com