

# BULLET ICE MAKER LTICE8-2



# BULLET ICE MAKER LTICE8-2

Used in Beverage industries, Biochemical and biological studies, Refrigerated perishables, Food industries, Medical laboratories, Healthcare, Pharmaceutical.

Also known as Ice Maker, Laboratory Bullet Ice Maker.

## LTICE8-2 BULLET ICE MAKER

Vertical stirring engine model

Refrigeration system produces uniform shape ice bullets

Bell-shaped ice with a hexagon hole

Highly durable stainless steel 304 2B built-up

Automatic controller increases working reliability

Large ice storage bin

Eco-friendly, energy saving non-Freon refrigerant



## SPECIFICATIONS

Model	LTICE8-2
Capacity	25 kg/24 h
Storage Bin Capacity	7 kg
Refrigerant	R134a
Condenser Type	Air Cooled
Overall Dimension	400x543x670 mm
Packing Dimension	470x615x730 mm
Weight	30/33 kg
Power	240 W
Power Supply	220V /110 V, 50 / 60 Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name
2900808008	Water filter



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)