

TURBIDITY METER LTMET10-1



TURBIDITY METER LTMET10-1

LTMET10-1 TURBIDITY METER



SPECIFICATIONS

Model	LTMET10-1
Measuring Unit	
Range	(0 ~20.00) NTU, (20.0 ~200.0) NTU
Resolution	0.01NTU, 0.1NTU, 1NTU
Accuracy	≤ ±8 %
Fluctuation	≤ ±1.0% FS/30 min
General	
Power Supply	AA Battery, AC Adapter 200-240 VAC input DC5V output
Dimension (mm) & Weight (kg)	100x220x80,0.8



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com