

TURBIDITY METER LTMET10-3



TURBIDITY METER LTMET10-3

LTMET10-3 TURBIDITY METER



SPECIFICATIONS

Model	LTMET10-3
Measuring Unit	
Range	(0.000 ~20.00) NTU, (20.01 ~200.0) NTU, (200 ~2000) NTU, (2001 ~4000) NTU
Accuracy	$\leq \pm 6\%$
Repeatability	$\leq 0.5 \%$
Fluctuation	$\leq \pm 0.5 \%$ FS/30xmin
General	
Power Supply	AC Adapter, 200-240 VAC input
Dimension (mm) & Weight (kg)	358x323x160, 8



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com