

# PH METER LTMET8-1



# PH METER LTMET8-1

## LTMET8-1 PH METER



### SPECIFICATIONS

Model	LTMET8-1
pH	
Range	-2.00 to 14.00 pH
Resolution	0.01 pH
Accuracy	± 0.01 pH
Calibration Points	Up to 2
Standard Recognition	NIST buffers
mV	
Range	-1800 to 1800 mV
Resolution	1 mV
Accuracy	± 0.1 %FS
Conductivity	
Range	/
Resolution	/
Accuracy	/
Reference Temperature	/
Measurement	
Reading Mode	Continuous
Reading Prompts	Reading, Stable
Temp. Compensation	MTC
Inputs	
pH Electrode	BNC(Q9)
Conductivity Probe	/
Display Options	
Backlight	Yes
Auto Shutdown	300, 600, 1200, 1800, 3600 sec., off
IP Rating	IP54
General	

Power	AC Adapter, 100-240 V AC input, DC 9V output
Dimensions	200x160x63 mm
Weight	600 g (1.32 lb)



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)