

PH METER LTMET8-5



PH METER LTMET8-5

LTMET8-5 PH METER



SPECIFICATIONS

Model	LTMET8-5
pH	
Range	-2.00 to 20.00 pH
Resolution	0.1, 0.01 pH
Accuracy	± 0.02 pH
Calibration Points	Up to 5
Standard Customization	Yes
mV	
Range	-2000.0 to 2000.0 mV
Resolution	0.1
Accuracy	± 0.3 mV or ± 0.1 % of reading whichever is greater
Temperature	
Range	- 5 to 110 °C, 23 to 230 °F
Unit	°C, °F
Resolution	0.1
Relative Accuracy	±0.2
Measurement	
Reading Mode	AutoRead(Fast, Medium, Slow), Timed, Continuous
Reading Prompts	Reading, Stable, Locked
Temp. Compensation	ATC, MTC
Data Management	
Data Storage	500 results each
GLP Features	Yes
Inputs	
pH Electrode	BNC(Q9)
Temp./DO Probe	4-pin aviation connector
Outputs	
USB	PC,printer

Display Options	
Backlight	Yes
Auto Shutdown	300, 600, 1200, 1800, 3600sec, off
IP Rating	IP65
Date and Time	Yes
General	
Power	Rechargeable Lithium batter; AC Adapte,100-240 V AC input, DC5V output
Dimensions	80 x 255 x 35 mm
Weight	400 g(0.88 lb)



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com