

# BIOLOGICAL MICROSCOPE LTMIC13-2



# BIOLOGICAL MICROSCOPE LTMIC13-2

Used in University, Medical, Laboratory, Research, Pharmaceutical.

## LTMIC13-2 BIOLOGICAL MICROSCOPE

Excellent optical function with infinite optical system, sientopf binocular viewing head, inclined at 30°, 360° and quadruple nosepiece. Convenient operation and special design for viewing incubating cell tissue. Valid mould proof design. It can be used under the high temperature and high humidity environment. Infinite optical system



## SPECIFICATIONS

Model	LTMIC13-2
Compensation Free Trinocular Head 30°	Standard
Viewing Head	
Siedentopf binocular viewing head, inclined at 30°, 360° rotatable, Interpupillary 48~75 mm	Standard
Siedentopf trinocular viewing head, inclined at 30°, 360° rotatable, Interpupillary 48~75 mm	Optional
Eyepiece	
WF 10 X/ 18	NA
P16 X/ 11	NA
WF10 X/ 20	Standard
WF15 X/ 16	Optional
Objective	
Achromatic objective 4X, 10X, 40X, 100X	NA
Plan achromatic objective 4X, 10X, 40X, 100X	NA
Infinite plan achromatic objective 20X, 60X	Optional
Infinite plan achromatic objective 4x, 10x, 40x, 100x	Standard
Nosepiece	
Quadruple nosepiece	NA
Backward quadruple nosepiece	Standard
Backward quintuple nosepiece	Optional
Focusing	
Coaxial coarse & fine adjustment, fine division 0.002 mm, coarse stroke 37.7 mm per rotation, fine stroke 0.2 mm per rotation, moving range 20 mm	Standard
Stage	
Double layers mechanical stage: 132x142 mm, moving range: 75x40 mm	NA
Rectangular 230x150 mm stage, moving range: 78x54 mm, using low positioned X/Y coaxial control knob	Standard
Condenser	
Abbe condenser NA1.25 with aperture diaphragm	Standard

Illumination	
S-LED Illumination, brightness adjustable	Standard
12 V/ 20 W Halogen lamp, brightness adjustable	Optional
Kohler illumination	Optional
Power Supply	
110~220 V, 50/60 Hz	Standard
Accessories	
Photo attachment, video attachment	Optional
Dark field condenser	NA
Polarization set, phase contrast kit	Optional
Package Size	
323x264x440 mm (MSC21-19), 380x260x510 mm (BMIC-602)	Standard
Gross Weight	
6.5 kg (MSC21-19), 10 kg (MSC21-20)	Standard



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)