

# BIOLOGICAL MICROSCOPE LTMIC14-4



# BIOLOGICAL MICROSCOPE LTMIC14-4

Used in University, Medical, Laboratory, Research, Pharmaceutical.

## LTMIC14-4 BIOLOGICAL MICROSCOPE

The integration of structure design to save space. The Microscope has a dual viewing head where one tube is 30° inclined, another tube is vertical and 360° rotatable. Valid mould proof design. It can be used under the high temperature and high humidity environment.



## SPECIFICATIONS

Model	LTMIC14-4
Viewing Head	Dual viewing head, one tube 30°, inclined and another tube vertical, 360° rotatable
Achromatic Objectives	4X, 10X, 40X (Spring) ,100X (Spring, Oil)
Eyepiece	Wide Field Eyepiece:WF10 X, (WF16 X optional)
Stage	Double Layer Mechanical Stage Size 120x125 mm
Focusing	Coaxial Coarse and Fine Adjustment, Focusing Range 30 mm, Focusing Interval 0.002 mm
Condenser	Abbe NA=1.25 with Iris Diaphragm & Filter
Illumination	LED(1 W) Lamp 220 V or 110 V, (LED 3 W lamp optional)
Power Supply	AC110/220 V ± 10 %, 50/60 Hz
Package Size	450x360x240 mm
Gross Weight	5 kg



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)