

STEREO ZOOM MICROSCOPE LTMIC16-1



STEREO ZOOM MICROSCOPE LTMIC16-1

Used in Pharmaceutical, University, Medical, Laboratory, Research.

LTMIC16-1 STEREO ZOOM MICROSCOPE

Excellent optical function with binocular viewing head inclined at 60°. Innovative stand structure and sharp image display. Convenient operation and special design for viewing incubating cell tissue. Valid mould proof design. It can be used under the high temperature and high humidity environment.



Product Image Coming Soon

SPECIFICATIONS

| Model | LTMIC16-1 |
|--|-----------|
| Viewing Head | |
| Binocular view head inclined at 60° | Yes |
| Trinocular view head inclined at 45° | No |
| Digital view head inclined at 45° | No |
| Eyepiece | |
| High point wide field adjustable eyepiece WF10x22 mm | Yes |
| Zoom Objective | |
| 0.67 X—4.5 X | Yes |
| Zoom Ratio | |
| 6.7:1 | Yes |
| Auxiliary Objective | |
| 0.5X /177 mm, 1.5X /47 mm, 2X/26 mm | Yes |
| Working Distance | |
| 115 mm | Yes |
| Focusing Bracket | |
| Focusing frame and adjustment knob | Yes |
| Interpupillary Distance | |
| 55-75 mm / 52-75 mm | Yes |
| Up Illumination | |
| LED lamp | Yes |
| Down Illumination | |
| LED lamp | Yes |
| Power Supply | |
| AC110/220 V ± 10%, 50/60 Hz | Yes |
| Package | |
| Cartoon with foam | Yes |

| | |
|--------------|----------------|
| Package Size | 500x260x370 mm |
| Gross Weight | 8 kg |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com