

STEREO ZOOM MICROSCOPE LTMIC16-2



STEREO ZOOM MICROSCOPE LTMIC16-2

Used in Pharmaceutical, University, Medical, Laboratory, Research.

LTMIC16-2 STEREO ZOOM MICROSCOPE

Excellent optical function with tinocular viewing head inclined at 60°. Innovative stand structure and sharp image display. Convenient operation and special design for viewing incubating cell tissue. Valid mould proof design. It can be used under the high temperature and high humidity environment. Trinocular view head inclined at 45°



SPECIFICATIONS

Model	LTMIC16-2
Viewing Head	
Binocular view head inclined at 60°	No
Trinocular view head inclined at 45°	Yes
Digital view head inclined at 45°	Yes
Eyepiece	
High point wide field adjustable eyepiece WF10x22 mm	Yes
Zoom Objective	
0.67 X — 4.5 X	Yes
Zoom Ratio	
6.7:1	Yes
Auxiliary Objective	
0.5X /177 mm, 1.5X /47 mm, 2X/26 mm	Yes
Working Distance	
115 mm	Yes
Focusing Bracket	
Focusing frame and adjustment knob	Yes
Interpupillary Distance	
55-75 mm / 52-75 mm	Yes
Up Illumination	
LED lamp	Yes
Down Illumination	
LED lamp	Yes
Power Supply	
AC110/220 V ± 10%, 50/60 Hz	Yes
Package	
Cartoon with foam	Yes

Package Size	500x260x370 mm
Gross Weight	8 kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com