

STEREO ZOOM MICROSCOPE LTMIC16-4



STEREO ZOOM MICROSCOPE LTMIC16-4

Used in Pharmaceutical, University, Medical, Laboratory, Research.

LTMIC16-4 STEREO ZOOM MICROSCOPE

Excellent optical function with Tranilluminator / inclined light illumination. Innovative stand structure and sharp image display. Convenient operation and special design for viewing incubating cell tissue. Valid mould proof design. It can be used under the high temperature and high humidity environment.



SPECIFICATIONS

Model	LTMIC16-4
Illumination	Tranilluminator / inclined light
View Head And Body	Inclined viewing Up-rotatable body
Objective	2X&4X, 90°rotated
Working Distance	57 mm
Eyepiece	WF10X/20
Power Supply	AC110/220 V ± 10%, 50/60 Hz
Package Size	360x170x290 mm
Gross Weight	6 kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com