

# STEREO ZOOM MICROSCOPE LTMIC16-5



# STEREO ZOOM MICROSCOPE LTMIC16-5

Used in Pharmaceutical, University, Medical, Laboratory, Research.

## LTMIC16-5 STEREO ZOOM MICROSCOPE

Excellent optical function with Natural light illumination. Innovative stand structure and sharp image display. Convenient operation and special design for viewing incubating cell tissue. Valid mould proof design. It can be used under the high temperature and high humidity environment.



### SPECIFICATIONS

Model	LTMIC16-5
Illumination	Natural light
View Head And Body	Inclined viewing Body 360°rotatable
Objective	2X&4X, 360°rotated
Working Distance	100 mm
Eyepiece	WF10X/20
Power Supply	AC110/220 V ± 10%, 50/60 Hz
Package Size	400x250x360 mm
Gross Weight	8 kg



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)