

MICROCIRCULATION MICROSCOPE LTMIC8-1



MICROCIRCULATION MICROSCOPE LTMIC8-1

Used in Microbiology Lab, Pharmaceutical, Hospitals, Medical, Universities, Scientific research, Institute.

Also known as Ultra microscope, High-magnification microscope, High resolution microscope.

LTMIC8-1 MICROCIRCULATION MICROSCOPE

Excellent optical function with infinite optical system.


Innovative stand structure and sharp image display.

Valid mould proof design, convenient operation and special design for viewing incubating cell tissue.

It can be used under the high temperature and high humidity environment.

Focus System: Coaxial coarse/ Slightly adjusted, Z axis movement, Mechanical upper and lower limits.

Illumination: LED illumination with high brightness and lower consumption, new luminaire with lower temperature and 10,000-200,000h working life, 90°revolve left and right.



Product Image Coming Soon

SPECIFICATIONS

Model	LTMIC8-1
Magnification	380X
Display Size	161x121 mm
Objective	10X
Stage	360°revolve, X,Y axis movement ,moving range 30x80 mm
Camera	Color camera system, HD LCD display
Package Size	270x340x380 mm
Gross Weight	5.5 kg
Power Supply	12 V



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com