

MICROCIRCULATION MICROSCOPE LTMIC8-2



MICROCIRCULATION MICROSCOPE LTMIC8-2

Used in Microbiology Lab, Pharmaceutical, Hospitals, Medical, Universities, Scientific research, Institute.

Also known as Ultra microscope, High-magnification microscope, High resolution microscope.

LTMIC8-2 MICROCIRCULATION MICROSCOPE

Excellent optical function with infinite optical system.

Innovative stand structure and sharp image display.

Valid mould proof design, convenient operation and special design for viewing incubating cell tissue.

It can be used under the high temperature and high humidity environment.

Focus System: Coaxial coarse/ Slightly adjusted, Z axis movement, Mechanical upper and lower limits.

Illumination: LED illumination with high brightness and lower consumption, new luminaire with lower temperature and 10,000-200,000h working life, 90°revolve left and right.



Product Image Coming Soon

SPECIFICATIONS

| | |
|---------------|---|
| Model | LTMIC8-2 |
| Magnification | 480X |
| Display Size | 204x153 mm |
| Objective | 10X |
| Stage | 360°revolve, X,Y axis movement ,moving range 30x80 mm |
| Camera | Color camera system, HD LCD display |
| Package Size | 300x380x380 mm |
| Gross Weight | 7 kg |
| Power Supply | 12 V |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com