

# BIOLOGICAL MICROSCOPE LTMIC9-4



# BIOLOGICAL MICROSCOPE LTMIC9-4

Used in University, Medical, Laboratory, Research, Pharmaceutical.

## LTMIC9-4 BIOLOGICAL MICROSCOPE

The Microscope has a Ball bearing quadruple nosepiece. It comes with an optional plan achromatic objective. It can be used under the high temperature and high humidity environment.



## SPECIFICATIONS

Model	LTMIC9-4
Viewing Head	
Sliding Trinocular Head Inclined at 45°	No
Compensation Free Trinocular Head inclined at 30°	No
10° LCD	Yes
Nosepiece	
Ball Bearing Quadruple Nosepiece	Yes
Eyepiece	
WF10X/ 18mm	Yes
Objective	
Achromatic 4X, 10X, 40X(S) 100X(S,O)	Yes
Plan Achromatic 4X, 10X, 40X(S) 100X(S,O)	Optional
Stage	
Double Layer Mechanical Stage 135x140mm	Yes
Condenser	
ABBE NA1.25 condenser with Iris Diaphragm & filter, rack&pinion adjustable	Yes
Illumination	
Build-in Illumination, Halogen lamp 6 V / 20 W, Power supply 90-230 V	Yes
Power Supply	
AC110/220 V +10 % , 50/60 Hz	Yes
Package	
Cartoon with foam	Yes
Package Size	
250x360x460 mm	Yes
Gross Weight	
8 kg	Yes



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)