

VORTEX MIXER LTMIX11-3



VORTEX MIXER LTMIX11-3

Used in Mixing, Laboratory, Research, Medical.

Also known as Laboratory Vortex Mixer.

LTMIX11-3 VORTEX MIXER

Touch operation or continuous mode

Used for various mixing applications with optional adapters

Various adapters fit for different tubes

Specially designed vacuum suction foot for body stability

Robust aluminum-cast construction



SPECIFICATIONS

Model	LTMIX11-3
Speed Range	25000 rpm
Permissible Ambient Temperature	5-40°C
Permissible Relative Humidity	80% RH
Orbital Diameter	4 mm
53	Motor rating output [W]
Operation	Continuous/ touch operation
Overall Dimension	127Wx130Hx160Dmm
Motor Type	Shaded pole motor
Protection Class	IP21
Weight	3.5 kg
Power	60 W
Power Supply	110-120 V/220-230 V, 50/60 Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
3301110008	Rubber vacuum suction foot with screw; 1pc	
3301110009	Platform pad	For φ99mm tubes and small vessels
3301110010	Standard top cup	For φ30mm tubes and small vessels
3301110011	Tube adapter	For 48 holes test tubes; $\varphi</math>6mm; used with Mx1027/Mx1026/PS1.1$

Accessory Code	Name	Description
3301110013	Tube adapter	For 12 holes test tubes; $\varphi</math>12mm ; used with Mx1027/Mx1026/PS1.1$
3301110014	Tube adapter	For 8 holes test tubes; $\varphi</math>16mm; used with Mx1027/Mx1026/PS1.1$
3301110015	Tube adapter	For 8 holes test tubes; $\varphi</math>20mm; used with Mx1027/Mx1026/PS1.1$

3301110012	Tube adapter	For 15 holes test tubes; φ10mm; used with Mx1027/ Mx1026/PS1.1
------------	--------------	----------------------------------------------------------------------

3301110016	Universal top plate	100mm plate attachment ; free of Rubber vacuum suction foot (18900158; 4pcs)
3301110017	Tube holding rod	Used with Mx1028 to Mx1032



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com