

# VORTEX MIXER LTMIX8-2



# VORTEX MIXER LTMIX8-2

Used in Mixing, Laboratory, Research, Medical.

Also known as Laboratory Vortex Mixer.

## LTMIX8-2 VORTEX MIXER

The max. diameter of the tubes is up to 30mm. Suit for mixing tubes /centrifugal tubes with prominent mixing effect

Outer shell and test tube base are made of plastic which is resistant to chemical attack. The thick base makes operation very steady and reliable

12V low voltage power adapter built in, smart and portable



## SPECIFICATIONS

|                         |                         |
|-------------------------|-------------------------|
| Model                   | LTMIX8-2                |
| Speed Range             | 4000 rpm                |
| Maximum Sample Capacity | 50 ml                   |
| Shaking Orbit           | 4.5 mm (circle)         |
| Switch mode             | Jog switch              |
| Power Input /Output     | 1.2 W / 0.8 W           |
| Voltage                 | DC12 V 1A               |
| Power                   | 12 W                    |
| Dimension               | W.110 x D.110 x H.86 mm |
| Net Weight              | 1.1 kgs                 |

# LABTARE

**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)