

# VORTEX MIXER LTMIX8-3



# VORTEX MIXER LTMIX8-3

Used in Mixing, Laboratory, Research, Medical.

Also known as Laboratory Vortex Mixer.

## LTMIX8-3 VORTEX MIXER

Suit for mixing tubes / centrifugal tubes with prominent mixing effect

Outer shell and touch concave are made of plastic which is resistance to chemical attack

The thick base makes operation very steady and reliable

12V low voltage power adapter built in, very smart and portable

Can adjust speed to meet different experimental requests



## SPECIFICATIONS

Model	LTMIX8-3
Speed Range	0-2500 rpm
Permissible Ambient Temperature	5-40°C
Permissible Relative Humidity	0.8
Orbital Diameter	4 mm
Motion	Orbital
58	Motor rating output [W&#93
Operation	Continuous/ touch operation
Overall Dimension	127Wx130Hx160D mm
Speed Display	Scale
Protection Class	IP21
Weight	3.5 kg
Power	60 W
Power Supply	110-120V /220-230V, 50/60Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
3301110008	Rubber vacuum suction foot with screw; 1pc	
3301110009	Platform pad	For <math>\varphi</math>99mm tubes and small vessels
3301110010	Standard top cup	For <math>\varphi</math>30mm tubes and small vessels
3301110011	Tube adapter	For 48 holes test tubes; $\varphi</math>6mm; used with Mx1027/ Mx1026/PS1.1$

Accessory Code	Name	Description
3301110013	Tube adapter	For 12 holes test tubes; $\varphi</math>12mm ; used with Mx1027/ Mx1026/PS1.1$
3301110014	Tube adapter	For 8 holes test tubes; $\varphi</math>16mm; used with Mx1027/ Mx1026/PS1.1$
3301110015	Tube adapter	For 8 holes test tubes; $\varphi</math>20mm; used with Mx1027/ Mx1026/PS1.1$

3301110012	Tube adapter	For 15 holes test tubes; φ10mm; used with Mx1027/ Mx1026/PS1.1
------------	--------------	--

3301110016	Universal top plate	100mm plate attachment ; free of Rubber vacuum suction foot (18900158; 4pcs)
3301110017	Tube holding rod	Used with Mx1028 to Mx1032



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)