

VORTEX MIXER LTMIX8-6



VORTEX MIXER LTMIX8-6

Used in Mixing, Laboratory, Research, Medical.

Also known as Laboratory Vortex Mixer.

LTMIX8-6 VORTEX MIXER

The perfect mixing radius combined with excellent 2D mixing technology and adjustable mixing speed ensures fast and perfect mixing of multiple samples.



SPECIFICATIONS

Model	LTMIX8-6
Speed Range	300 ~3000 rpm
Touch Vortex Speed	3500 rpm
Shaking Orbit	3 mm
Time Range	1 s~99 h 59 min
Working Noise	<50dB
Voltage	AC 100~230V, 50/60Hz
Power	45 W
Fuse	250 V, 1 A, ϕ 5 x 20
Dimension	W.196 x D.270 x H.170 mm
Net Weight	7.7 kgs

ACCESSORIES

Accessory Code	Name	Capacity
3300813008	Tube Stand	96x0.2ml tube stand. Applied to 0.2ml PCR tubes; PCR strips and 96-wellPCR plates(semi-skirted and unskirted)
3300813009	Tube Stand	4x0.5ml Tube stand
3300813010	Tube Stand	24x1.5/2.0ml Tube stand



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com