

# MINI MIXER LTMIX9-1



# MINI MIXER LTMIX9-1

Used in Mixing, Laboratory, Research, Medical.

Also known as Laboratory Vortex Mixer.

## LTMIX9-1 MINI MIXER

Mini and compact appearance, powerful function, sucker type under-chassis, super shockproof, ideal for high speed operation

Multiple mixing modes, touch-operation, continuous operation and timing set

Speed adjustable, and wide speed range, Stable low speed, powerful high speed

Shaking blocks convenient for replacement, stable and reliable, sturdy and durable eccentric bearing designed

Various tube holders are for optional (Suitable for Eppendorf Tubes)



## SPECIFICATIONS

Model	LTMIX9-1
Speed Range	0 ~ 2500 rpm
Shaking Orbit	4 mm
Shaking Mode	Circle
Operation	knob + scale
Operation Mode	continuous or press running
Voltage	AC 220V / AC 110V, 50/60 Hz
Power	60 W
Fuse	250 V, 1 A, $\Phi$ 5x20
Dimension	W.132 x D.156 x H.150 mm
Net Weight	2.5 kgs

## ACCESSORIES

Accessory Code	Name	Description
3301008022	Standard head	Suit for tubes and containers with dia. < 30mm
3301008023	Standard platform	Standard platform to hold foam tube holder

## OPTIONAL ACCESSORIES

Accessory Code	Name
3301008008	Convex pad
3301008009	Foam tube holder
3301008010	Foam tube holder
3301008011	Foam tube holder
3301008012	Foam tube holder

Accessory Code	Name
3301008015	Tube holder platform
3301008016	Tube holder platform
3301008017	Plate holder platform
3301008018	Universal platform for tube stands and all standard plates

3301008013	Foam tube holder
3301008014	Tube holder platform

3301008019	Universal platform for tube stands and all standard plates
3301008020	Universal platform for tube stands and all standard plates
3301008021	Customized



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)