

4°C BLOOD BANK REFRIGERATOR LTREF10-5



4°C BLOOD BANK REFRIGERATOR LTREF10-5

Used in Immunohematology Department, Blood Bank, Laboratory, Research, Medical, Pharmaceutical.

Also known as Laboratory Blood Bank Refrigerator.

LTREF10-5 4°C BLOOD BANK REFRIGERATOR

450ml blood bag x 200



SPECIFICATIONS

Model	LTREF10-5
Capacity	368 L
Temperature Range	+4°C±1
Ambient Temperature	10~32°C
Bloodbags(450ml)	200
Refrigerant	R134a
Defrost Method	Auto Defrost
Controller	Microprocessor Control with digital Display
Control Sensor	NTC
Cooling Type	Forced air cooling
Compressor brand/QTY	EMBRACO/1
Shelf/drawer	Shelves/4
Alarms	High/low temperature,Power off, Temperature controller malfunction, Backup battery malfunction, 8h Power failure backup system(alarm)
Insulation	PURF
Material Inside	Stainless steel
Material Outside	Color sprayed steel
Internal Size	685×429×1380 mm
External Size	785×565×1920 mm
Package Size	890×650×2100 mm
CBM	1.21
20GP/40GP/40HQ	21/44/44
NT./GT. (kg)	128/135
Power consumption (kW.h/24h)	4.28
Power	325W
Current (A)	2.2A

Voltage/frequency (V/Hz)	220/50Hz
Standard Accessories	Caster, 1 Heating & Foam glass door type outdoor, 4 Inner door, Test hole 25mm, Thermal Printer, LED Light,
Optional Accessory	Temperature recorder Datalogger



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com