

PORTABLE FREEZER LTREF15-2



PORTABLE FREEZER LTREF15-2

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

LTREF15-2 PORTABLE FREEZER



SPECIFICATIONS

Model	LTREF15-2
External Size	335x225x205 mm
Internal Size	295x167x168 mm
Inner Capacity	8 L
Product Configuration	1pc refrigerator, 4 pcs Crystal storage ice box 400ml, 1pc digital thermometer
Insulation Material	EPS
Gross Weight	1.8 kg
Packing Size	380x260x290 mm



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com