

PORTABLE FREEZER LTREF15-3



PORTABLE FREEZER LTREF15-3

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

LTREF15-3 PORTABLE FREEZER



SPECIFICATIONS

Model	LTREF15-3
External Size	425x262x305 mm
Internal Size	325x170x225 mm
Inner Capacity	12 L
Product Configuration	1pc refrigerator, 1pc tray 1pc partition 6 pcs Crystal storage ice box 400ml 1pc digital thermometer
Insulation Material	PU
Gross Weight	4 kg
Packing Size	450x330x330 mm



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com