

PORTABLE FREEZER LTREF15-5



PORTABLE FREEZER LTREF15-5

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

LTREF15-5 PORTABLE FREEZER



SPECIFICATIONS

| | |
|-----------------------|--|
| Model | LTREF15-5 |
| External Size | 450x330x300 mm |
| Internal Size | 370x235x240 mm |
| Inner Capacity | 20 L |
| Product Configuration | 1pc refrigerator, 1pc strap 1pc tray 9 pcs Crystal storage ice box 400ml 1pc digital thermometer 1pc partition |
| Insulation Material | PU |
| Gross Weight | 5 kg |
| Packing Size | 490x350x330 mm |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com