

PORTABLE FREEZER LTREF15-6



PORTABLE FREEZER LTREF15-6

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

LTREF15-6 PORTABLE FREEZER



SPECIFICATIONS

| | |
|-----------------------|--|
| Model | LTREF15-6 |
| External Size | 530x325x340 mm |
| Internal Size | 480x270x265 mm |
| Inner Capacity | 33 L |
| Product Configuration | 1pc refrigerator, 12 pcs Crystal storage ice box 400ml 1pc digital thermometer |
| Insulation Material | PU |
| Gross Weight | 5.5 kg |
| Packing Size | 550x340x360 mm |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com