

# PORTABLE REFRIGERATOR LTREF16-1



# PORTABLE REFRIGERATOR LTREF16-1

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

## LTREF16-1 PORTABLE REFRIGERATOR

- 1.The storage container can keep content cooling up to 24 hours
- 2.Insulation material is PU CFC-free polyurethane
- 3.Made of environmental PP material,food grade, non-toxic
- 4.Storage container is durable for vaccine storage, blood products storage,dry ice storage,frozen specimen storag etc



## SPECIFICATIONS

Model	LTREF16-1
External Size (WxDxH)	340x230x280 mm
Internal Size (WxDxH)	226x140x210 mm
Capacity	6 L
Temperature Range	2°C-8°C
Storage Time (Ambient 30°C)	24 h
Material	
External	PP
Internal	PU
Ice Vessels	6(CMB-400)
Max Loading	6 kg
Optional Accessory	LED display thermometer
Net Weight	1.9 kg
Gross Weight	5 kg
Package Size (WxDxH)	360x250x300 mm



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)