

# PORTABLE REFRIGERATOR LTREF16-3



# PORTABLE REFRIGERATOR LTREF16-3

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

## LTREF16-3 PORTABLE REFRIGERATOR

- 1.The storage container can keep content cooling up to 24 hours
- 2.Insulation material is PU CFC-free polyurethane
- 3.Made of environmental PP material,food grade, non-toxic
- 4.Storage container is durable for vaccine storage, blood products storage,dry ice storage,frozen specimen storag etc



## SPECIFICATIONS

Model	LTREF16-3
External Size (WxDxH)	440x300x340 mm
Internal Size (WxDxH)	350x220x250 mm
Capacity	17 L
Temperature Range	2°C-8°C
Storage Time (Ambient 30°C)	24 h
Material	
External	PP
Internal	PU
Ice Vessels	2(CMB-1200)
Max Loading	15 kg
Optional Accessory	LED display thermometer
Net Weight	3.2 kg
Gross Weight	6.5 kg
Package Size (WxDxH)	460x330x370 mm



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)