

PORTABLE REFRIGERATOR LTREF16-4



PORTABLE REFRIGERATOR LTREF16-4

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

LTREF16-4 PORTABLE REFRIGERATOR

- 1.The storage container can keep content cooling up to 24 hours
- 2.Insulation material is PU CFC-free polyurethane
- 3.Made of environmental PP material,food grade, non-toxic
- 4.Storage container is durable for vaccine storage, blood products storage,dry ice storage,frozen specimen storag etc



SPECIFICATIONS

Model	LTREF16-4
External Size (WxDxH)	440x270x310 mm
Internal Size (WxDxH)	375x215x250 mm
Capacity	20 L
Temperature Range	2°C-8°C
Storage Time (Ambient 30°C)	18 h
Material	
External	PP
Internal	PU
Ice Vessels	8(CMB-400)
Max Loading	18 kg
Optional Accessory	LED display thermometer
Net Weight	3.5 kg
Gross Weight	7.5 kg
Package Size (WxDxH)	460x290x330 mm



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com