

PORTABLE REFRIGERATOR LTREF16-5



PORTABLE REFRIGERATOR LTREF16-5

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

LTREF16-5 PORTABLE REFRIGERATOR

- 1.The storage container can keep content cooling up to 24 hours
- 2.Insulation material is PU CFC-free polyurethane
- 3.Made of environmental PP material,food grade, non-toxic
- 4.Storage container is durable for vaccine storage, blood products storage,dry ice storage,frozen specimen storag etc



SPECIFICATIONS

Model	LTREF16-5
External Size (WxDxH)	570x310x370 mm
Internal Size (WxDxH)	450x225x300 mm
Capacity	33 L
Temperature Range	2°C-8°C
Storage Time (Ambient 30°C)	20 h
Material	
External	PP
Internal	PU
Ice Vessels	3(CMB-1200)
Max Loading	10 kg
Optional Accessory	LED display thermometer
Net Weight	4.5 kg
Gross Weight	9.5 kg
Package Size (WxDxH)	590x330x390 mm



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com