

PORTABLE REFRIGERATOR LTREF16-6



PORTABLE REFRIGERATOR LTREF16-6

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

LTREF16-6 PORTABLE REFRIGERATOR

- 1.The storage container can keep content cooling up to 24 hours
- 2.Insulation material is PU CFC-free polyurethane
- 3.Made of environmental PP material,food grade, non-toxic
- 4.Storage container is durable for vaccine storage, blood products storage,dry ice storage,frozen specimen storag etc



SPECIFICATIONS

Model	LTREF16-6
External Size (WxDxH)	615x395x450 mm
Internal Size (WxDxH)	525x300x355 mm
Capacity	54 L
Temperature Range	2°C-8°C
Storage Time (Ambient 30°C)	24 h
Material	
External	PP
Internal	PU
Ice Vessels	4(CMB-1500)
Max Loading	11 kg
Optional Accessory	LED display thermometer
Net Weight	7.8 kg
Gross Weight	17 kg
Package Size (WxDxH)	800x550x510 mm



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com