

# 4°C BLOOD BANK REFRIGERATOR LTREF8-4



# 4°C BLOOD BANK REFRIGERATOR LTREF8-4

Used in Immunohematology Department, Blood Bank, Laboratory, Research, Medical, Pharmaceutical.  
Also known as Laboratory Blood Bank Refrigerator.

## LTREF8-4 4°C BLOOD BANK REFRIGERATOR

Vertical cabinet

Microprocessor-based temperature controller, Temperature range  $4 \pm 1^\circ\text{C}$

Large screen LCD display of temperature, display accuracy of  $\pm 0.1^\circ\text{C}$

Polyurethane insulation is 80-95 mm

Forced air circulation system, defrost automatic

Digital display of upper and lower temperature

Audible and visual alarm

Exterior made from high quality steel plate, Interior is stainless steel

3-layer transparent door with lock, door with heater around to prevent ice leak.

5 Inner drawer (stainless steel)

Interior fluorescent light

The high-effective air-cooled condenser and the fin-type evaporator

Standard: Temperature printer

Optional: Chart recorder



## SPECIFICATIONS

Model	LTREF8-4
Capacity	358 L
Temperature Range	$+4^\circ\text{C}$ , $+2 \sim +10^\circ\text{C}$
Refrigerant	R134a, CFC free
Refrigeration System	International famous compressor and Germany EBM fan motor.
Defrost Method	No frost design
Controller	Microprocessor Controller with LCD Display
Control Sensor	2
Alarms	High & low temperature alarm, No battery alarm, Door ajar alarm, Filter blocking alarm, System failure alarm
Shelves	5 units of drawers
External Size	700x628x1962 mm
Package Size	512x515x1270 mm
Power Supply	220V,50/ 60Hz, 110V,50/60Hz, 220V,50Hz/ 120V,60Hz

# LABTARE

**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)