

# +4/+22 BLOOD TRANSPORT BOX LTREF8-7



# +4/+22 BLOOD TRANSPORT BOX LTREF8-7

Used in Immunohematology Department, Blood Bank, Laboratory, Research, Medical, Pharmaceutical.  
Also known as Laboratory Blood Bank Refrigerator.

## LTREF8-7 +4/+22 BLOOD TRANSPORT BOX

4°C for storing blood

22°C for storing blood platelet

Audible and visual alarm

Four units of baskets inside as standard, which can load 50 blood bags of 450ml



## SPECIFICATIONS

Model	LTREF8-7
Capacity	90 L
Temperature Range	4°C or 22°C
Ambient Temperature	-5°C or 43°C
Compressor	Danfoss compressor
Alarms	High & low temperature alarm
Insulation	50 mm
Material	Whole LLDPE
External Size	1000Lx565Wx862H mm
Package Size	1135Lx636Wx1069H mm
Optional Accessories	Hydraulic support bar, Remote alarm system
Weight	73/90 kg
Power Supply	DC12/24V/ AC100-240V 50/60Hz Automatic selection



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)