

# TUBE ROTATOR LTROT10-2



# TUBE ROTATOR LTROT10-2

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Rotator, Laboratory Rotator.

## LTROT10-2 TUBE ROTATOR

The rotator comes with a disk type design

Digital control speed and allowing procedures to be accurate

LCD displays timer

End-over-end or rolling action, gives gentle but effective mixing

Adjustable speed range of 10 to 70rpm

Adjustable mixing angle between 0 and 90 °

Used with a wide variety of accessories

Continuous and timed quiet operation



## SPECIFICATIONS

Model	LTROT10-2
Speed Range	10-70 rpm
Permissible Ambient Temperature	5-40°C
Permissible Relative Humidity	0.8
Tilt Angle (°)	0-90
Operation	Continuous/ timed operation
Timer	Yes
Timing Range	1-1199 min
Dimension	510Dx190Wx260H mm
Motor Type	DC motor
Protection Class	IP21
Display	LCD
Weight	8 kg
Power	40 W
Power Supply	100-240V, 50/60Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
3501009008	Disk accessory	For 1.5ml x 60 centrifuge tubes holder
3501009009	Disk accessory	For 15ml x 16 centrifuge tubes holder
3501009010	Disk accessory	For 50ml x 8 centrifuge tubes holder
3501009011	Disk support rods	4 pcs; use with BX-103920; BX-103921; BX-103922



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)