

TUBE ROTATOR LTROT10-4



TUBE ROTATOR LTROT10-4

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Rotator, Laboratory Rotator.

LTROT10-4 TUBE ROTATOR

The rotator comes with a rotisserie type design

Digital control speed and allowing procedures to be accurate

LCD displays timer

End-over-end or rolling action, gives gentle but effective mixing

Adjustable speed range of 10 to 70 rpm

Used with a wide variety of accessories

Continuous and timed quiet operation



SPECIFICATIONS

Model	LTROT10-4
Speed Range	10-70 rpm
Permissible Relative Humidity	0.8
Tilt Angle (°)	90
Operation	Continuous/ timed operation
Timer	Yes
Timing Range	1-1199 min
Dimension	510Dx190Wx260H mm
Motor Type	DC motor
Protection Class	IP21
Display	LCD
Weight	8 kg
Power	40 W
Power Supply	100-240V, 50/60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
3501011008	Rotisserie accessory	For 1.5ml x 48 centrifuge tubes held horizontally
3501011009	Rotisserie accessory	For 15ml x 24 centrifuge tubes held horizontally
3501011010	Rotisserie accessory	For 50ml x 24 centrifuge tubes held horizontally
3501011011	Rotisserie accessory	For 1.5ml x 32 centrifuge tubes held vertically
3501011012	Rotisserie accessory	For 15ml x 16 centrifuge tubes held vertically
3501011013	Rotisserie accessory	For 50ml x 16 centrifuge tubes held vertically



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com