

3D ROTATOR LTROT11-1



3D ROTATOR LTROT11-1

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Rotator, Laboratory Rotator.

LTROT11-1 3D ROTATOR

Unique combination of adjustable rotating frame.

A variety of rotation modes, from mild type to somersault type rotation mode.

Positive and negative rotation just gently press the rotating frame when instrument is working.

Widely used in molecular biology, histo-chemistry, bio-chemistry, clinical applications, etc.

Provide three types of rotating frame.



SPECIFICATIONS

| | |
|-------------|--------------------------|
| Model | LTROT11-1 |
| Speed Range | 18 rpm |
| Voltage | AC 220V/AC 110V, 50/60Hz |
| Power | 4 W |
| Dimension | 240x145x210 mm |
| Net Weight | 1 kg |

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com