

# VERTICAL ROTATOR LTROT12-2



# VERTICAL ROTATOR LTROT12-2

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Rotator, Laboratory Rotator.

## LTROT12-2 VERTICAL ROTATOR

Smart appearance, solid structure, simple operation

LED display time and speed

DC brushless motor, stable performance, low noise and long using life

360° vertical circulation rotation, operating smooth and steady

Different tube holder frames for optional to suit different types of tubes



## SPECIFICATIONS

Model	LTROT12-2
Speed Range	10 ~ 80 rpm
Shaking Orbit	circle
Time Range	1min~99h59min
Voltage	AC 100~230V, 50/60Hz
Power	16 W
Fuse	250V, 1A, $\Phi$ 5x20
Dimension	W.470 x D.190 x H.160mm
Net Weight	2.8 kgs

## ACCESSORIES

Accessory Code	Name	Capacity
3501209008	Tube Holder frame	1.5mlx64
3501209009	Tube Holder frame	15mlx22
3501209010	Tube Holder frame	50mlx16



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)