

# VERTICAL ROTATOR LTROT12-3



# VERTICAL ROTATOR LTROT12-3

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Rotator, Laboratory Rotator.

## LTROT12-3 VERTICAL ROTATOR

Humanized designed operation keyboard

LED display time and speed

Several optional turntable supplied, suit for the tubes from 1.5ml to 50ml

The turntable's gradient can be adjusted from 0 to 90 degree

Timing function

Beep alarm when program time ends

DC brushless motor, stable performance, low noise and long using life



## SPECIFICATIONS

Model	LTROT12-3
Speed Range	10 ~ 80 rpm
Gradient Angle Range	0°~90°
Time Range	1min~9999min
Voltage	AC 100~230V, 50/60 Hz
Power	35 W
Fuse	250V, 1A, $\Phi$ 5x20
Dimension	W.206 x D.260 x H.300mm
Net Weight	4.5 kgs

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Capacity
3501210008	Turntable for rotary mixing solution	50ml V-bottom tube x 10	
3501210009	Turntable for rotary mixing solution	1.5ml tube x 64	
3501210010	Turntable for rotary mixing solution	15ml V-bottom tube x 18	
3501210011	Foam tube holder		50ml centrifuge tubex12 15ml centrifuge tubex12 1.5ml centrifuge tubex30
3501210012	Turntable		50ml centrifuge tubex10 15ml centrifuge tubex10 1.5ml centrifuge tubex30
3501210013	Double turntable connector		



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)