

# TUBE ROTATOR LTROT8-2



# TUBE ROTATOR LTROT8-2

Used in Mixing, Laboratory, Research, Medical.

Also known as Rotator, Laboratory Rotator.

## LTROT8-2 TUBE ROTATOR

It can be used in the humid environment such as refrigerator and refrigerated cabinet.

Several optional turntable supplied, suit for the tubes from 1.5ml to 50ml.

One button to start/stop, easy to operate, can run for a long time, running smoothly without noise at 40 rpm.



## SPECIFICATIONS

Model	LTROT8-2
Speed Range	40 rpm
Gradient Angle Range	45° (fixed angle)
Voltage	AC 220V/AC 110V, 50/60 Hz
Power	14 W
Fuse	250V, 1A, $\Phi$ 5x20
Dimension	W.260 x D.300 x H.233mm
Net Weight	1.8 kgs

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Capacity
3501210008	Turntable for rotary mixing solution	50ml V-bottom tube x 10	
3501210009	Turntable for rotary mixing solution	1.5ml tube x 64	
3501210010	Turntable for rotary mixing solution	15ml V-bottom tube x 18	
3501210011	Foam tube holder		50ml centrifuge tubex12 15ml centrifuge tubex12 1.5ml centrifuge tubex30
3501210012	Turntable		50ml centrifuge tubex10 15ml centrifuge tubex10 1.5ml centrifuge tubex30
3501210013	Double turntable connector		



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)