

# TUBE ROLLER LTROT9-2



# TUBE ROLLER LTROT9-2

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Roller, Laboratory Roller.

## LTROT9-2 TUBE ROLLER

Rocking and rolling motions for complete mixing

LCD display the speed and time

Speed range of 10-70 rpm

Six roller design with small footprint and easy clean

Designed for continuous or timed operation



## SPECIFICATIONS

Model	LTROT9-2
Speed Range	10-70 rpm
Roller Length	279 mm
Amplitude	24 mm
Maximum Load Capacity	4 kg
Operation	Continuous or timed
Dimension	444Lx279Wx125H mm
Motor Type	DC brushless
Weight	5.5 kg
Power Supply	100-220V, 50/60Hz



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)