

# TUBE ROLLER LTROT9-3



# TUBE ROLLER LTROT9-3

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Roller, Laboratory Roller.

## LTROT9-3 TUBE ROLLER

Compact and firm designed, simple operation, LED display time and speed.

Two types selectable - rolling type and rolling & rolling type and rolling & swinging type. Continuous working mode, safe and stable.

With stable performance motor, operation speed adjustable, max. speed to 80rpm.

Low noise, long using life.

Can compatible with different size of tubes. Large range of application.



## SPECIFICATIONS

Model	LTROT9-3
Speed Range	5 ~ 80 rpm
Operation Mode	Rolling and Fluctuation
Amplitude	25 mm
Roller Quantity	6
Roller Length	280 mm
Max. Load	4 kgs
Time Range	1min~99h59min
Voltage	AC 100~230V, 50/60 Hz
Power	25 W
Fuse	250V, 1A, $\Phi$ 5x20
Dimension	W.450 x D.260 x H.120 mm



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)