

TUBE ROLLER LTROT9-4



TUBE ROLLER LTROT9-4

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Roller, Laboratory Roller.

LTROT9-4 TUBE ROLLER

Compact and firm designed, simple operation, LED display time and speed.

Two types selectable - rolling type and rolling & swinging type. Continuous working mode, safe and stable.

With stable performance motor, operation speed adjustable, max. speed to 80rpm.

Low noise, long using life.

Can compatible with different size of tubes. Large range of application.



SPECIFICATIONS

| | |
|-----------------|--------------------------|
| Model | LTROT9-4 |
| Speed Range | 5 ~ 80 rpm |
| Operation Mode | Rolling |
| Roller Quantity | 6 |
| Roller Length | 280 mm |
| Max. Load | 4 kgs |
| Time Range | 1min~99h59min |
| Voltage | AC 100~230V, 50/60 Hz |
| Power | 25 W |
| Fuse | 250V, 1A, Φ 5x20 |
| Dimension | W.450 x D.260 x H.120 mm |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com