

DRYING OVEN OVE12-040



DRYING OVEN OVE12-040

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

OVE12-040 DRYING OVEN

Digital display all 304 stainless steel inner chamber and outer

Intelligent PID control system

Large observation window

Tight buckle handle



SPECIFICATIONS

Model	OVE12-040
Capacity	40 L
Temperature Accuracy	±1°C
Temperature Fluctuation	±1%
Temperature Uniformity	0.5°C
Temperature Resolution	0.1°C
Time for Heating	≤6 min (20°C to 57°C)
Timer Range	0-999 min
Internal Dimension (WxDxH)	350×350×350 mm
Overall Dimension (WxDxH)	640×550×530 mm
Weight	35.5/38.5 kg
Power	1000 W
Power Supply	220 V, 50 Hz



Labtare USA

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@labtare.com | Website: labtare.com