

# DRYING OVEN OVE12-070



# DRYING OVEN OVE12-070

Heating Oven, Laboratory Drying Oven

Used in Drying, Roasting, Laboratory, Research, Engineering, Industry, Agricultural, Biology.

Also known as Laboratory Oven, Laboratory general Purpose Oven.

## OVE12-070 DRYING OVEN

Micro-computer intelligent control system with LED display

Forced air convection maintains temperature uniformity

High quality stainless steel inner chamber with electrostatic spray treated exterior

Large observation window

Easy to clean semicircular corners with adjustable shelves

Silicon rubber sealing ring prevents heat loss

Over temperature alarm with timer



## SPECIFICATIONS

|                            |                 |
|----------------------------|-----------------|
| Model                      | OVE12-070       |
| Capacity                   | 70 L            |
| Temperature Range          | RT+5-250°C      |
| Temperature Accuracy       | ±1°C            |
| Temperature Uniformity     | ±4.5°C at 250°C |
| Temperature Resolution     | ±1°C            |
| Timer Range                | 0-999 min       |
| Internal Dimension (WxDxH) | 420x450x350 mm  |
| Overall Dimension (WxDxH)  | 630x740x550 mm  |
| Package Size               | 750x820x630 mm  |
| Shelves                    | 2               |
| Weight                     | 44.5/47.5 kg    |
| Power                      | 1200 W          |
| Power Supply               | 220V, 60Hz      |

## OPTIONAL ACCESSORIES

| Accessory Code | Name  | Description |
|----------------|---|-------------|
| LS54389        | Intelligent program temperature controller  |             |
| LS54401        | Independent limit of temperature controller |             |





**Labtare USA**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)