

DRYING OVEN OVE12-070T



DRYING OVEN OVE12-070T

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

OVE12-070T DRYING OVEN

Digital display all 304 stainless steel inner chamber and outer

Intelligent PID control system

Large observation window

Tight buckle handle



SPECIFICATIONS

Model	OVE12-070T
Capacity	70 L
Temperature Accuracy	$\pm 1^{\circ}\text{C}$
Temperature Fluctuation	$\pm 1\%$
Time for Heating	≤ 6 min (20°C to 57°C)
Timer Range	0-999 min
Temperature control accuracy	$\leq \pm 0.5^{\circ}\text{C}$
Internal Dimension (WxDxH)	450X350X450 mm
Overall Dimension (WxDxH)	740x550x630 mm
Display Accuracy	0.1°C
Weight	45/48 kg
Power	1200 W
Power Supply	220 V, 50 Hz



Labtare USA

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@labtare.com | Website: labtare.com