

DRYING OVEN OVE12-130



DRYING OVEN OVE12-130

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

OVE12-130 DRYING OVEN

Digital display all 304 stainless steel inner chamber and outer

Intelligent PID control system

Large observation window

Tight buckle handle



SPECIFICATIONS

Model	OVE12-130
Capacity	130 L
Temperature Accuracy	±1°C
Temperature Fluctuation	±1%
Timer Range	0-999 min
Internal Dimension (WxDxH)	550x450x550 mm
Overall Dimension (WxDxH)	890x630x750 mm
Weight	60/64 kg
Power	1600 W
Power Supply	220 V, 50 Hz



Labtare USA

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@labtare.com | Website: labtare.com