

# MAGNETIC STIRRER STR17-010



# MAGNETIC STIRRER STR17-010

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

## STR17-010 MAGNETIC STIRRER

Speed range of 0-1500 rpm

Maximum stirring quantity 10 L

Glass ceramic plate, excellent chemical resistance and easy to clean



## SPECIFICATIONS

|             |            |
|-------------|------------|
| Model       | STR17-010  |
| Capacity    | 10 L       |
| Speed Range | 0-1500 rpm |

## OPTIONAL ACCESSORIES

| Accessory Code | Name                | Description                       | Dimension  |
|----------------|---------------------|-----------------------------------|------------|
| LS55170        | Stirrer bar         |                                   | 10mm x 6mm |
| LS55171        | Stirrer bar         |                                   | 15mm x 8mm |
| LS55172        | Stirrer bar         |                                   | 20mm x 8mm |
| LS55173        | Stirrer bar         |                                   | 25mm x 8mm |
| LS55174        | Stirrer bar         |                                   | 30mm x 6mm |
| LS55175        | Stirrer bar         |                                   | 40mm x 8mm |
| LS55176        | Stirrer bar         |                                   | 50mm x 8mm |
| LS55177        | Stirrer bar remover | Length of 200mm                   |            |
| LS55195        | Stirrer bar         | Stirrer bar (65mm x 8mm), 1pc/pk  |            |
| LS55196        | Stirrer bar         | Stirrer bar (80mm x 13mm), 1pc/pk |            |





**Labtare USA**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)