

# MAGNETIC HOTPLATE STIRRER STR21-340T



# MAGNETIC HOTPLATE STIRRER STR21-340T

Rugged and compact magnetic hotplate stirrers are ergonomically designed and provide uniform heating surfaces. It delivers reliable and consistent stirring quickly and efficiently. The low-profile design of these magnetic hotplate stirrers makes them easy to store and fit into fume hoods.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Hotplate Stirrer, Digital Magnetic Stirrer.

## STR21-340T MAGNETIC HOTPLATE STIRRER

Available for heating and stirring standard/non-standard reaction flasks between 50 ml and 20 L

DC brushless motor make stable operation, low noise, long using life and no spark

Metal plate shell, firm and durable, thermostability and anti-corrosion

The control panel designed 30° slopes for convenient operation

Steady working at low speed, powerful working at high speed

Temperature Sensor:PT1000



## SPECIFICATIONS

Model	STR21-340T
Temperature Setting Range	30°C~340°C
Temperature Stability	±3°C
1	Speed Range
RT~340°C	Permissible Ambient Temperature
5°C	Temperature Increment
350°C	Control Accuracy with Sensor PT1000
50°C	Timing Range
1 min to 99h 59 min	Motor rating Input [W&#93
Enamel	Platform Dimension
Φ137 mm	Speed Display
156Wx270Dx70H mm	Fuse
250 V,4 A/8 A, Φ5x20	Weight
2.2 kg	Power
600 W	Power Supply

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS55169	Support rod and Temperature probe	





**Labtare USA**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)