

BULLET ICE MAKERS



BULLET ICE MAKERS

Used in Beverage industries, Biochemical and biological studies, Refrigerated perishables, Food industries, Medical laboratories, Healthcare, Pharmaceutical.

Also known as Ice Maker, Laboratory Bullet Ice Maker.

LTICE8 BULLET ICE MAKER

- Vertical stirring engine model
- Refrigeration system produces uniform shape ice bullets
- Bell-shaped ice with a hexagon hole
- Highly durable stainless steel 304 2B built-up
- Automatic controller increases working reliability
- Large ice storage bin
- Eco-friendly, energy saving non-Freon refrigerant



SPECIFICATIONS

Model	LTICE8-1	LTICE8-2	LTICE8-3	LTICE8-4
Capacity	20 kg/24 h	25 kg/24 h	50 kg/24 h	80 kg/24 h
Storage Bin Capacity	6 kg	7 kg	14 kg	25 kg
Refrigerant	R134a			
Condenser Type	Air Cooled			
Overall Dimension	330x470x605 mm	400x543x670 mm	400x510x825 mm	500x612x900 mm
Packing Dimension	405x555x665 mm	470x615x730 mm	470x590x875 mm	560x670x950 mm
Weight	28/32 kg	30/33 kg	35/38 kg	51/56 kg
Power	220 W	240 W	380 W	580 W
Power Supply	220V /110 V, 50 / 60 Hz			

OPTIONAL ACCESSORIES

Accessory Code	Name
2900808008	Water filter



LTICE8-1



LTICE8-2



LTICE8-3



LTICE8-4

LTICE10-12 BULLET ICE MAKER

- Transparent Bullet Ice
- Air Cooling
- Stainless Steel 304 Exterior
- Adjustable Feet
- Slide Away Access Door
- World Famous Brand Compressor
- Fully Automatic Controller
- Improved Design For Higher Reliability



SPECIFICATIONS

Model	LTICE10-12
Capacity	120 kg/24 h
Storage Bin Capacity	55 kg
Refrigerant	R134a
Condenser Type	Air Cooled
Overall Dimension	680x611x1050 mm
Packing Dimension	750x675x1150 mm
Weight	70/75 kg
Power	800 W
Power Supply	220 V /110 V, 50 / 60 Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com