

CUBE ICE MAKERS



CUBE ICE MAKERS

Used in Ice Cuber Machine, Laboratory Cube Ice Maker.

LTICE10 CUBE ICE MAKER

High grade stainless steel SUS304 is used which is anticorrosive and durable. The compact structure is simple and space-saving.

Fully automatically computer controlled, convenient for use.

Automatic detection and protection in case of operation, water deficiency, ice full and malfunction.

Highly efficient fluoride-free compressor, energy saving and environment-friendly.

Suitable for Pubs, KTV and dancing halls, gyms, cold drinks shops, milk tea shops, cafes, convenience stores, restaurants and other entertainment venues.

Split-type structure, large ice bin, large ice storage capacity and strong applicability.



SPECIFICATIONS

Model	LTICE10-13	LTICE10-14	LTICE10-15	LTICE10-16
Capacity	200 kg/24 h	300 kg/24 h	500 kg/24 h	1000 kg/24 h
Storage Bin Capacity	145 kg	210 kg	280 kg	205 kg
Refrigerant	R134a / R404a			
Condenser Type	Water cooled			-
Machine Head	565x600x560 mm	760x600x560 mm	760x600x740 mm	1230x630x740 mm
Cabinet	565x800x900 mm	760x800x980 mm		1230x920x1095 mm
Packing Machine Head	650x640x850 mm	845x705x850 mm	845x705x950 mm	1320x680x950 mm
Packing Cabinet	600x850x980 mm	800x850x1060 mm		1290x1000x1160 mm
Weight	71/79 kg	98/107 kg	123/133 kg	205/225 kg
Power Supply	220 V, 60 Hz			380 V, 60 Hz



LTICE10-13



LTICE10-14



LTICE10-15



LTICE10-16

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com