

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINETS



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINETS

It is specially designed storage unit which helps to segregate and store weak acids and alkalis. It is constructed from high-quality materials like corrosion resistant Stainless steel or polyethylene which prevents chemical reactions that might arise within acids or alkalis. The cabinet specifically has separate compartments or shelves with designated areas for acids and alkalis to prevent cross-contamination and ensure the safe storage practices.

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research.

Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

LTCAB8 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



SPECIFICATIONS

Model	LTCAB8-1	LTCAB8-2	LTCAB8-3	LTCAB8-4
Capacity	4 gallon (15L)	12 gallon (45L)	22 gallon (83L)	30 gallon (114L)
Material	Double layer anti fire cold-rolled steel PP tray			
Adjustable Shelf	1 pc		3 pcs	1 pc
Door	Single / Manual			Double / Manual
External Size	430x430x560 mm	590x460x890 mm	600x460x1650 mm	1090x460x1120 mm
Package Size	550x560x725 mm	710x580x1066 mm	720x580x1815 mm	1210x680x1285 mm
Gross Weight	38 Kg	66 Kg	116 Kg	117 Kg

Model	LTCAB8-5	LTCAB8-6	LTCAB8-7
Capacity	45 gallon (170L)	60 gallon (227L)	90 gallon (340L)
Material	Double layer anti fire cold-rolled steel PP tray		
Adjustable Shelf	2 pcs		
Door	Double / Manual		
External Size	1090x460x1650mm	860x860x1650 mm	1090x860x1650 mm
Package Size	1210x580x1815 mm	980x980x1815 mm	1210x980x1815 mm
Gross Weight	155 Kg	195 Kg	243 Kg



LTCAB8-1



LTCAB8-2



LTCAB8-3



LTCAB8-4



LTCAB8-5



LTCAB8-6



LTCAB8-7

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com